



Part **Basket without Lock**
 Part Number **02830-03**
 Serial Number **P11 001001 – P11 001117**



Test Certificate

		USP / Ph.Eur. Tolerances		
Dimensions	Basket	Basket rim O.D. Ø	25.0 mm	± 3.0 mm
		Height	37.0 mm	± 3.0 mm
		Wire Ø	0.25 mm	0.22 – 0.31mm
		Wire openings	0.4 mm	0.36 - 0.44 mm
		Screen opening height	27.0 mm	± 1.0 mm
		Screen O.D. Ø	22.2 mm	± 1.0 mm
		Clear opening Ø	20.2 mm	± 1.0 mm
Material	Stainless Steel (electrolytic polished)			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopoeia requirements (at date of approval, see below).

Navi Mumbai, 2022-07-18

18-07-2022 15:33:49

Moreshwar Pulekar
 SOTAX Quality Management
 C-64 | First Floor | Suchi Business Park | TTC Industrial area | MIDC Turbhe | Navi Mumbai 400705 | India