



Part Basket without Cover, 40 Mesh  
 Part Number 02830  
 Serial Number N11020283 - N11020354



# Test Certificate

		USP / Ph.Eur. Tolerances		
<b>Dimensions</b>	Basket	Basket Lower Ring O.D. Ø	25.0 mm	± 3.0 mm
		Height	37.0 mm	± 3.0 mm
		Wire Ø	0.25 mm	0.22 – 0.31mm
		Wire openings	0.4 mm	0.36 - 0.44 mm
		Screen opening height	27.0 mm	± 1.0 mm
		Screen O.D. Ø	22.2 mm	± 1.0 mm
		Clear opening Ø	20.2 mm	± 1.0 mm
<b>Material</b>	Stainless Steel (electrolytic polished)			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements of a random sampling of the serial number range are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Westborough, 2024-03-21



21.03.2024 13:07:24

Jeff Schiavoni  
 SOTAX Quality Management  
 SOTAX Corp | 2400 Computer Dr | Westborough, MA 01581 | USA