



Part Paddle with Stirring Shaft
 Part Number 07525-04
 Serial Number N16010652 - N16010777



Test Certificate

| | | | | USP / Ph.Eur. Tolerances |
|-------------------|---|-------------------|---------------|-----------------------------|
| Dimensions | Paddle | Height | 19.0 mm | ± 0.5 mm |
| | | Upper edge length | 74.5 mm | ± 0.5 mm |
| | | Lower edge length | 42.0 mm | ± 1.0 mm |
| | | Thickness | 4.0 mm | ± 1.0 mm |
| | | Radius | 1.2 mm | ± 0.2 mm |
| | | Radius | 41.5 mm | ± 1.0 mm |
| | Shaft | Diameter | 9.4 – 10.1 mm | -- |
| Material | Stainless Steel (electrolytic polished) | | | |

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Aesch, 2022-11-23

Bruno Balli
 SOTAX Quality Management
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland