


Test Certificate

Part: **Glass Container inactinic brown**

Number: 09126-2

Serial Number: 72853 - 72857

For SOTAX CTS Content Testing Systems



Dimensions:	Inside \varnothing	145 mm	± 2 mm
	Height	164.5 mm	± 2 mm

With straight bottom

Material: Borosilicate glass

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed.

SOTAX Quality Management
Bruno Balli
2017-11-20



Signature