


Test Certificate

Part: **Vessel Mini Stirrer**
 Number: 10583
 Serial Number: 71840 - 71874

For SOTAX AT Dissolution Testing Systems



Dimensions:				
A	Outside Ø	124.5 mm	+ 0.3 / - 0.2 mm	
B	Inside Ø	100.0 mm	± 0.5 mm	
C	Outside Ø	50.0 mm	± 1.0 mm	
D	Inside Ø	40.0 mm	± 1.0 mm	
E	Height	211.5 mm	± 0.5 mm	

With hemispherical bottom

Material: Borosilicate glass

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed.

SOTAX Quality Management
 Bruno Balli
 2017-08-21


 Signature