


Test Certificate

Part: **Nozzle 25 mm**
Number: 231126

For SOTAX PF1



| | | | |
|-------------|-----------------|----------|---------------------------|
| Dimensions: | Upper Inside Ø | 30.2 mm | USP/Ph.Eur. Tolerances |
| | Bottom Inside Ø | 25.00 mm | 0 / - 0.1 mm |
| | Chamfer | 1.0 mm | ± 0.01 mm |
| | Funnel (Cone) | 40 ° | ± 0.2 mm |
| | | | ± 0.8 ° |

Material: Stainless Steel (electrolytic polished)

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at time of delivery).

SOTAX Quality Management
Bruno Balli
2018-08-13


Signature