



Part Paddle with Suction Opening  
 Part Number 02985-2  
 Serial Number N13004525 – N13004625



# Test Certificate

|                   |   |                   |               | USP / Ph.Eur.<br>Tolerances |
|-------------------|---|-------------------|---------------|-----------------------------|
| <b>Dimensions</b> | Paddle                                  | Height            | 19.0 mm       | ± 0.5 mm                    |
|                   |   | Upper edge length | 74.5 mm       | ± 0.5 mm                    |
|                   |   | Lower edge length | 42.0 mm       | ± 1.0 mm                    |
|                   |   | Thickness         | 4.0 mm        | ± 1.0 mm                    |
|                   |   | Radius            | 1.2 mm        | ± 0.2 mm                    |
|                   |   | Radius            | 41.5 mm       | ± 1.0 mm                    |
|                   | Shaft                                   | Diameter          | 9.4 – 10.1 mm | --                          |
| <b>Material</b>   | Stainless Steel (electrolytic polished) |                   |               |                             |

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Aesch, 2025-01-20

Buner Cornelius  
 SOTAX Quality Management  
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland