



Part Paddle with Stirring Shaft
 Part Number 07525-04
 Serial Number N16011473 - N16011623



Test Certificate

			USP / Ph.Eur. Tolerances	
Dimensions	Paddle	Height	19.0 mm	± 0.5 mm
		Upper edge length	74.5 mm	± 0.5 mm
		Lower edge length	42.0 mm	± 1.0 mm
		Thickness	4.0 mm	± 1.0 mm
		Radius	1.2 mm	± 0.2 mm
		Radius	41.5 mm	± 1.0 mm
	Shaft	Diameter	9.75 mm	± 0.35 mm
Material	Stainless Steel (electrolytic polished)			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Aesch, 2023-03-03

Tayyib Özudogru
 SOTAX Quality Management
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland