




Part 1 Liter Vessel inactinic brown
Part Number 14537-02
Serial Number 98207 - 98213



| | | | Tolerances | |
|------------|--------------------|---------------------------|-------------------|----------|
| Dimensions | Vessel | Inside Ø a | 102.9 to 105.1 mm | -- |
| | | Height b | 186.9 mm | ± 0.7 mm |
| | | With hemispherical bottom | | |
| Material | Borosilicate glass | | | |

Aesch, 2024-07-12

Tayyib Özüdogru
SOTAX Quality Management
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland

Test Certificate